



Bioreactors/Fermenters

	ambr 15 cell culture	ambr 250 modular	ambr 250 high throughput	Biostat B	Biostat B-DCU	Biostat Cultibag RMI1	Biostat STR
Microbial fermentations	No	Yes	Yes	Yes	Yes	Yes	Yes
Mammalian fermentations	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Working volume (minimum)	15 ml (10 ml)	250 ml (100 ml)	250 ml (100 ml)	1L (0.4L)	1L (0.4L)	1L (0.1-0.5L)	50L (12.5L)
				2L (0.6L)	2L (0.6L)	2L (0.2-1.0L)	200L (50L)
				5L (0.6L)	5L (0.6L)	10L (1-5L)	500L (125L)
				10L (1.5L)	10L (1.5L)	20L (2-10L)	1000L (250L)
						50L (5-25L)	2000L (500L)
						100L (10-50L)	
						200L (20-100L)	
Max. number of vessels	24 or 48	12 or 24	2-8	2	6	2	2
Twin-option	N/A	N/A	N/A	Yes	No	Yes	Yes
Single use	Yes	Yes	Yes	Possible	Yes	Yes	Yes